



# Important Info For Scanning



**DynaFlex® Accepts Scans  
From Any Scanner**

When scanning your patient, always keep in mind the desired appliance, visualizing areas that will be incorporated in its fabrication and use. This will ensure a proper turnaround time and accuracy in the fabricated appliance. Please refer to the chart below for our scanning requirements:

Appliance	Upper Palate Needed	Bite Needed	Include 2nd Molars	Margin Extensions	Special Instructions
*All Upper Fixed & Removable Appliances	YES	YES, for all functional appliances.	YES, If appliance design includes them.	7-10mm extension beyond the gingival margin.	Remove arch-wires if patient is in orthodontic braces.
All Lower Fixed & Removable Appliances	N/A	YES, for all functional appliances.	YES, If appliance design includes them.	7-10mm extension beyond the gingival margin.	Remove arch-wires if patient is in orthodontic braces.
EZ-Align®, Defend®, Essix	NO	NO	YES, If appliance design includes them.	3-4mm extension beyond the gingival margin.	Remove arch-wires if patient is in orthodontic braces.
Sleep Apnea Device	NO	YES	YES, If appliance design includes them.	7-10mm extension beyond the gingival margin.	N/A

**\*Excluding Bonded Retainer, Herbst or Mara which do not require major connectors.**

020518 © 2018 DynaFlex®, St. Louis, MO 63074. Printed in U.S.A. All rights reserved.

It is a violation of copyright law to reproduce all or part of this material, including photography, without the permission of DynaFlex®.