

Patient Cooperation Not Required



Class II Solution



The CS-2000[®] is Your Class II and Class III Solution

This revolutionary system will guarantee results with your Class II and Class III orthodontic cases. The CS-2000[®] Niti springs are a DynaFlex[®] exclusive, utilizing a patented, low-continuous force, closed coil spring and a specially built "key-hole" end which fits over a patented slotted base pivot similar to a Herbst[®] housing nut.



Class III Solution

070308 © Copyright 2008. Herbst[®] is a registered trademark of Dentaaurum. CS-2000[®] patent #6719557. Pivot with slotted base patent #6036488.



The slotted pivot allows for an arch wire to slide through the base, thereby reducing the overall profile of the attachment. The slotted pivot can be welded to a band or crown and also used as an arch wire attachment for amazing versatility.

A. CS-2000[®] Niti Spring (available in 10mm and 7mm) • B. Pivot with Slotted Base (.022 x .025) • C. Screw for Pivot • D. Wrench • E. Bur (optional)

All items sold in kit or individually

Since
1965



1-800-489-4020
www.dynaflex.com

DynaFlex[®]
Your Total Orthodontic Solution!SM