

# Crystal clear!

## Stain resistant!

Manufactured of monocrystalline sapphire, these brackets are crystal clear! ClearViz+™ brackets are a low profile bracket with a high bond strength and are stain resistant. Choose from Roth and MBT\* prescriptions.

## ClearViz+™

### Bracket System




Color-Coded  
Marking System

Made Of Pure  
Monocrystalline  
Sapphire

Almost  
Invisible  
Against  
Most Tooth  
Shades

Heat  
Polished  
For  
Smooth  
Edges

Base Is Able To Flex During  
Debonding, Providing  
More Patient Comfort



 **DynaFlex®**

web: [www.dynaflex.com](http://www.dynaflex.com) ph: 800-489-4020

041415 © 2014 DynaFlex®, St. Louis, MO 63074. Printed in U.S.A. All rights reserved. It is a violation of copyright law to reproduce all or part of this material, including photography, without the permission of DynaFlex®.