

Krystal™ Brackets are the choice for a **true ceramic bracket**

Manufactured from 99.9% pure alumina oxide, these polycrystalline brackets are non-staining and do not discolor. Krystal Brackets are a true twin design offering outstanding bond, high strength, precise tip and torque and are very easy to ligate. Your patients will be pleased with the comfort and aesthetic look!

Krystal™ Ceramic Bracket System

Impervious
To Staining

Made Of 99.9% Alumina Oxide

Low Profile Design

Ultra Smooth Slot For
Less Friction

Compound
Contour Base

Mechanical Lock Base For
Ease In Debonding



 **DynaFlex®**

web: www.dynaflex.com ph: 866-346-5665

062217 © 2017 DynaFlex®, St. Louis, MO 63074. Printed in U.S.A. All rights reserved. It is a violation of copyright law to reproduce all or part of this material, including photography, without the permission of DynaFlex®.