


The Revolution™ Self-Ligating System

In creating our latest innovation, the DynaFlex® engineers thought outside the box to produce a truly unique, passive, self-ligating bracket system that performs as smoothly as it looks. With a sleek, ultra low profile design for ultimate patient comfort, the Revolution™ bracket system is perhaps the easiest self-ligating system available on the market, requiring no special instruments or training of any kind. Simply bond it, clip it, and watch it work. It is an effective self-ligating bracket system that allows you to focus on the treatment, and not the appliances. It truly is revolutionary.



Classic Twin Design

The Revolution™ bracket system allows for fast and accurate bracket placement due to its familiar twin design. Even though Revolution™ is a self-ligating bracket system, the under tie wing area allows for the option of engaging elastomeric ligatures or power chain over the clip or after the clip is removed. Revolution™ is designed for easy removal of the self-ligating clip to enhance the finishing process. Available in Roth and MBT® prescriptions.

 **DynaFlex®**

web: www.dynaflex.com ph: 800-489-4020

102814 © 2014 DynaFlex®, St. Louis, MO 63074. Printed in U.S.A. All rights reserved. It is a violation of copyright law to reproduce all or part of this material, including photography, without the permission of DynaFlex®.