





























Atlas Mini Bracket System, Roth Rx - 10 Pack

MAXILLARY	REF	HK	TRQ	ANG	M/D	IN/OUT	.018 RIGHT	.018 LEFT	.022 RIGHT	.022 LEFT
Central	UR  UL 	-	12	5	3.4mm	0.61mm	#201-R8-11	#201-R8-21	#201-R2-11	#201-R2-21
Lateral	UR  UL 	-	8	9	2.8mm	0.80mm	#201-R8-12	#201-R8-22	#201-R2-12	#201-R2-22
Lateral w/ Hk	UR  UL 	D	8	9	2.8mm	0.80mm	#201-R8-12D	#201-R8-22D	#201-R2-12D	#201-R2-22D
Cuspid	UR  UL 	-	-2	11	3.0mm	0.58mm	#201-R8-13	#201-R8-23	#201-R2-13	#201-R2-23
Cuspid w/ Hk	UR  UL 	D	-2	11	3.0mm	0.58mm	#201-R8-131	#201-R8-231	#201-R2-131	#201-R2-231
1st Bicuspid	UR  UL 	-	-7	0	2.8mm	0.59mm	#201-R8-14	#201-R8-24	#201-R2-14	#201-R2-24
1st Bicuspid w/ Hk	UR  UL 	M	-7	0	2.8mm	0.59mm	#201-R8-141	#201-R8-241	#201-R2-141	#201-R2-241
1st Bicuspid w/ DG Hk	UR  UL 	D	-7	0	2.8mm	0.59mm	#201-R8-141D	#201-R8-241D	#201-R2-141D	#201-R2-241D
2nd Bicuspid	UR  UL 	-	-7	0	2.8mm	0.59mm	#201-R8-15	#201-R8-25	#201-R2-15	#201-R2-25
2nd Bicuspid w/ Hk	UR  UL 	M	-7	0	2.8mm	0.59mm	#201-R8-151	#201-R8-251	#201-R2-151	#201-R2-251
2nd Bicuspid w/ DG Hk	UR  UL 	D	-7	0	2.8mm	0.59mm	#201-R8-151D	#201-R8-251D	#201-R2-151D	#201-R2-251D
MANDIBULAR	REF	HK	TRQ	ANG	M/D	IN/OUT	.018 RIGHT	.018 LEFT	.022 RIGHT	.022 LEFT
Lower Anteriors	LR  LL 	-	0	0	2.3mm	0.78mm	#201-R8-31	Universal	#201-R2-31	Universal
Lower Anteriors Wide	LR  LL 	-	0	0	2.9mm	0.78mm	#201-R8-31S	Universal	#201-R2-31S	Universal
Lateral w/ Hk	LR  LL 	D	0	0	2.3mm	0.78mm	#201-R8-32D	#201-R8-42D	#201-R2-32D	#201-R2-42D
Cuspid	LR  LL 	-	-11	5	2.8mm	0.61mm	#201-R8-33	#201-R8-43	#201-R2-33	#201-R2-43
Cuspid w/ Hk	LR  LL 	D	-11	5	2.8mm	0.61mm	#201-R8-331	#201-R8-431	#201-R2-331	#201-R2-431
1st Bicuspid	LR  LL 	-	-17	0	2.8mm	0.64mm	#201-R8-34	#201-R8-44	#201-R2-34	#201-R2-44
1st Bicuspid w/ Hk	LR  LL 	M	-17	0	2.8mm	0.64mm	#201-R8-341	#201-R8-441	#201-R2-341	#201-R2-441
1st Bicuspid w/ DG Hk	LR  LL 	D	-17	0	2.8mm	0.64mm	#201-R8-341D	#201-R8-441D	#201-R2-341D	#201-R2-441D
2nd Bicuspid	LR  LL 	-	-22	0	2.8mm	0.69mm	#201-R8-35	#201-R8-45	#201-R2-35	#201-R2-45
2nd Bicuspid w/ Hk	LR  LL 	M	-22	0	2.8mm	0.69mm	#201-R8-351	#201-R8-451	#201-R2-351	#201-R2-451
2nd Bicuspid w/ DG Hk	LR  LL 	D	-22	0	2.8mm	0.69mm	#201-R8-351D	#201-R8-451D	#201-R2-351D	#201-R2-451D

Atlas Mini Bracket System, Roth Rx - 5x5 Set

FULL SETS	.018 SLOT	.022 SLOT
Upper/Lower 5x5 No Hks	#201-R8-NHK	#201-R2-NHK
Upper/Lower 5x5 Hks on 3's	#201-R8-3	#201-R2-3
Upper/Lower 5x5 Hks on 3,4,5's	#201-R8-345	#201-R2-345
Upper/Lower 5x5 Hks on 3,4,5's Distal Gingival Bicuspid	#201-R8-345DG	#201-R2-345DG

Atlas Mini Bracket System With Offset Base, Roth Rx - 10 Pack

MAXILLARY	REF	HK	TRQ	ANG	M/D	IN/OUT	.018 RIGHT	.018 LEFT	.022 RIGHT	.022 LEFT
1st Bicuspid	UR ●	-	-7	0	2.8mm	0.59mm	#201-R8-140-2	#201-R8-240-2	#201-R2-140-2	#201-R2-240-2
	UL ●									
1st Bicuspid w/ Hk	UR ●	D	-7	0	2.8mm	0.59mm	#201-R8-1410-2	#201-R8-2410-2	#201-R2-1410-2	#201-R2-2410-2
	UL ●									
2nd Bicuspid	UR ●	-	-7	0	2.8mm	0.59mm	#201-R8-150-2	#201-R8-250-2	#201-R2-150-2	#201-R2-250-2
	UL ●									
2nd Bicuspid w/ Hk	UR ●	D	-7	0	2.8mm	0.59mm	#201-R8-1510-2	#201-R8-2510-2	#201-R2-1510-2	#201-R2-2510-2
	UL ●									
MANDIBULAR	REF	HK	TRQ	ANG	M/D	IN/OUT	.018 RIGHT	.018 LEFT	.022 RIGHT	.022 LEFT
1st Bicuspid	LR ●	-	-17	0	2.8mm	0.64mm	#201-R8-340-2	#201-R8-440-2	#201-R2-340-2	#201-R2-440-2
	LL ●									
1st Bicuspid w/ Hk	LR ●	D	-17	0	2.8mm	0.64mm	#201-R8-3410-2	#201-R8-4410-2	#201-R2-3410-2	#201-R2-4410-2
	LL ●									
2nd Bicuspid	LR ●	-	-22	0	2.8mm	0.69mm	#201-R8-350-2	#201-R8-450-2	#201-R2-350-2	#201-R2-450-2
	LL ●									
2nd Bicuspid w/ Hk	LR ●	D	-22	0	2.8mm	0.69mm	#201-R8-3510-2	#201-R8-4510-2	#201-R2-3510-2	#201-R2-4510-2
	LL ●									

Atlas Mini Bracket System With Oval Base, Roth™ Rx - 10 Pack

MAXILLARY	REF	HK	TRQ	ANG	M/D	IN/OUT	.018 RIGHT	.018 LEFT	.022 RIGHT	.022 LEFT
1st Bicuspid	UR ●	-	-7	0	2.8mm	0.59mm	#201-R8OV-14	#201-R8OV-24	#201-R2OV-14	#201-R2OV-24
	UL ●									
1st Bicuspid w/ Hook	UR ●	D	-7	0	2.8mm	0.59mm	#201-R8OV-141	#201-R8OV-241	#201-R2OV-141	#201-R2OV-241
	UL ●									
2nd Bicuspid	UR ●	-	-7	0	2.8mm	0.59mm	#201-R8OV-15	#201-R8OV-25	#201-R2OV-15	#201-R2OV-25
	UL ●									
2nd Bicuspid w/ Hook	UR ●	D	-7	0	2.8mm	0.59mm	#201-R8OV-151	#201-R8OV-251	#201-R2OV-151	#201-R2OV-251
	UL ●									
MANDIBULAR	REF	HK	TRQ	ANG	M/D	IN/OUT	.018 RIGHT	.018 LEFT	.022 RIGHT	.022 LEFT
1st Bicuspid	LR ●	-	-17	0	2.8mm	0.64mm	#201-R8OV-34	#201-R8OV-44	#201-R2OV-34	#201-R2OV-44
	LL ●									
1st Bicuspid w/ Hook	LR ●	D	-17	0	2.8mm	0.64mm	#201-R8OV-341	#201-R8OV-441	#201-R2OV-341	#201-R2OV-441
	LL ●									
2nd Bicuspid	LR ●	-	-22	0	2.8mm	0.69mm	#201-R8OV35	#201-R8OV-45	#201-R2OV-35	#201-R2OV-45
	LL ●									
2nd Bicuspid w/ Hook	LR ●	D	-22	0	2.8mm	0.69mm	#201-R8OV351	#201-R8OV-451	#201-R2OV-351	#201-R2OV-451
	LL ●									